

Titrate Installation Guide

Table of Contents

<i>Document Control</i>	<i>1</i>
Revision History	1
Document Approval.....	1
<i>Objective</i>	<i>2</i>
<i>Prerequisites</i>	<i>2</i>
<i>Database Setup</i>	<i>3</i>
Installation	3
Create Database.....	5
<i>GIT Pull</i>	<i>6</i>
API	6
UI.....	6
<i>Run API</i>	<i>6</i>
Steps	6
<i>Run UI</i>	<i>8</i>
Steps	8
<i>Run Python Server</i>	<i>8</i>
Ubuntu Installation	9
Python Installation.....	14
Steps to run python code	15

Document Control

Revision History

Version	Date	Updated By	Description
V1.0	11-02-2024	Nikhil	Ubuntu & Python Installation
V1.1	08-03-2024	Ramesh	Database API & UI Setup

Document Approval

Date	Name	Role	Status
	Ramachandra Damera	Manager	Approved

Objective

This document provides a step-by-step guide to install and configure the application on a client environment. It covers the backend API and frontend UI setup, including database configuration using PostgreSQL and environment setup for both components.

➤ Prerequisites

- Git installed on the machine
- A 64-bit version of Windows 11
- VS Code
- Node.js and npm
- Ubuntu above 22.0
- Python 3.10+ and pip in Ubuntu
- PostgreSQL 13+ installed and running
- Access to Git repositories for API and UI
- An active internet connection

➤ Database Setup

Step 1: Connect to Your Server

- Launch **pgAdmin** and expand the **Servers** panel.
- Right-click your server (e.g., *PostgreSQL 14* or *localhost*) and ensure you're logged in as a role with **CREATE** permission (usually *postgres*).

Step 2: Open "Create Database" Dialog

- Under your server node, right-click **Databases** → select **Create** → **Database**.

Step 3: General Tab – Define Basics

- **Database:** Type your desired name (e.g., *my_first_database*).
- **Owner:** Choose the role (default: *postgres*).
- **Check the Port**
- **Give the host name, password and user name.**

Note: Remember the credentials given here.

Step 4: Save & Finish

- Click **Save** to create the database. The new entry will appear under the **Databases** node.

- Get the pull from GIT for API and UI respectively.

<https://github.com/Brase-Technologies-Pty-Ltd/Mitte-Analytics.git>

➤ Configuring and Running the API

- Open the file manager and navigate to *Titrate\api*
- Open Command Prompt from the navigated path.
- Open .env file and configure the data base credentials as noted while creating a database.



```
1 # Environment
2 NODE_ENV=development
3
4 # Database Configuration
5 PORT=4000
6 DB_USERNAME=postgres
7 DB_PASSWORD=postgres
8 DB_NAME=imprestStock_db
9 DB_HOST=localhost
10
11 # Real SMTP (production)
12 SMTP_HOST=mail.kloners.in
13 SMTP_PORT=465
14 SMTP_SECURE=true
15 SMTP_FROM=stock_isms@kloners.in
16 SMTP_USER=stock_isms@kloners.in
17 SMTP_PASS=Rameshaa@16
18
19 # Dev SMTP (Papercut - runs on localhost)
20 DEV_SMTP_HOST=13.126.10.72
21 DEV_SMTP_PORT=1025
22 DEV_SMTP_SECURE=false
23 DEV_SMTP_FROM=local@devmail.test
24
25 # Test customer email
26 CUSTOMER_EMAIL=ayyalarameshnaidu@gmail.com
27
28 # Default Admin Credentials
29 DEFAULT_ADMIN_EMAIL=admin@admin.com
30 DEFAULT_ADMIN_PASSWORD = Admin@123
31
32 # Purchase Order Folder
33 PO_FOLDER=
   "D:\Projects\Projects\brasetech\source\api\purchase_orders"
```

- Give the path for PO_FOLDER where you wish to save the file.
- Open the terminal and Run the command “npm install”.
- Run the command “*npm start*” in the command prompt. // command to start the API.

Note: The Titrate application will be running on the default port localhost: 4000

On running the API, default user will be created and added to the database. You can use the same to login and add new users and any other data respectively to the application.


DEFAULT_ADMIN_EMAIL=admin@admin.com

DEFAULT_ADMIN_PASSWORD = Admin@123

For confirmation you will be displayed with the message with port and IP address.

➤ Setting Up and Running the Web Server

- Go to file manager and navigate to *Titrate\api_ui*
- Open command prompt from the navigated path
- Open .env file and configure URL and App name and save.

A screenshot of a code editor with a dark background. It shows three lines of code in a .env file. Line 1: VITE_APP_API_BASE_URL. Line 2: = http://localhost:4000/api/. Line 3: VITE_APP_NAME = "Titrate". A blue bracket on the right side groups lines 1 and 2 together. The code is color-coded: VITE_APP_API_BASE_URL is blue, = is black, http://localhost:4000/api/ is black, VITE_APP_NAME is blue, = is black, and "Titrate" is orange.

```
1 VITE_APP_API_BASE_URL
  = http://localhost:4000/api/
2
3 VITE_APP_NAME = "Titrate"
```


- Run the command “npm install”.
- Run the command “*npm run dev*” in the command prompt // command to start the web

Note: The Titrate application will be running on the default port localhost: 80


➤ Run Python Server

Steps to run Python Code:


- Open Ubuntu and type “explorer.exe .”





- Pull the code from git repo
- Go to python folder




- Open the folder in VS python code you will be able to see the files on the left side.




- Now open extensions using command “Ctrl+Shift+x” and download WSL, python




- Click Green button at the left bottom of the screen to open a remote window





- A search filter will be opened at the top of screen. Select “Connect to WSL Distro....” and then select “ubuntu”







- It may take some time to set up VS Code in Ubuntu. Once completed, you will notice 'WSL: Ubuntu' displayed in the bottom-left corner of the window.




- Open and Select the folder and click “OK”




- Make sure that the all files are visible now on your left side navigation bar.



- Open config.ini which is present in files named “file_aus” and “file_austest”



- Change the Serial number and API key that matches with the camera running (If any).
- Change the paths of [Paths] in both the config files of file_aus and file_austest folders. For model_path copy and paste the path of best.pt from weights->weights->best.pt. And follow the same for all paths in [Paths].




```

1 [Camera]
2   serial_number = Q4EH-UMTG-L6LB
3   api_key = 82cfcc5f0967127091bcd855835f39aebaf4df97
4
5 [Paths]
6   model_path = /home/techrammy/mittedelivery/MitteAnalytics/runs/detect/train4/weights/best.pt
7   predict_path = /home/techrammy/mittedelivery/MitteAnalytics/predictpath
8   predict1_path = /home/techrammy/mittedelivery/MitteAnalytics/predictpath/predict
9   barcode_reader_exe = /home/techrammy/mittedelivery/MitteAnalytics/BarcodeReaderCLI/bin/BarcodeReaderCLI
10
11 [API]
12   get_api_url = http://13.126.223.104:2900/api/imprestproduct/fetch
13   get_api_url = http://13.126.223.104:2900/api/imprest/all
14   put_api_url = http://13.126.223.104:2900/api/imprestproduct/
15   post_api_url = http://13.126.223.104:2900/api/motionalert/objectDetected
16
17 [AWS]
18   access_key = MITTAVOUNCIVUFACTER

```


Note: Here, backend code is getting integrated with frontend using [API] section from config.ini file as in above image



- To open terminal select the ...-> Terminal->New Terminal from menu bar




NOTE: There are 2 ways to run the code,
a) Create a venv (virtual environment) & run initially for installing all required packages and then run the code

- Create a Virtual Environment







- Give permissions for the folder we have placed in by running command:
“chmod –R 777 (path/to/folder)”



- Run the command in the terminal “./final_au.sh” for running the camera of hospital impress and run command in new terminal “./final_austest.sh” for running camera of AU test impress.





```
root@amesh-ayyala:/home/techramy/mittedelivery/MitteAnalytics# ./final_austest.sh
WARNING: user config directory '/root/.config/Ultralytics' is not writeable, defaulting to '/tmp' or GID. Alternatively you can define a YOLO_CONFIG_DIR environment variable for this path.
AU TEST IMPREST

0: 384x640 2 shelves, 5 trays, 20 barcodes, 52 mboxes, 1 capsule, 152.9ms
Speed: 5.4ms preprocess, 152.9ms inference, 1.4ms postprocess per image at shape (1, 3, 384, 640)
Combined Results:
[]
Count not found for product ALL03
Count not found for product ALL04
Count not found for product AMIT01
Count not found for product AMIT1
Count not found for product AMIT2
```

This will get us the count and also will update the count in the UI Dashboard.